



EpiLife® Medium

Catalog Number: M-EPI-500-CA
500 ml

Instructions for storage and use

Product Description

EpiLife is a sterile, liquid tissue culture medium prepared with 60 µM calcium chloride for the growth of normal human epidermal keratinocytes and normal human corneal epithelial cells. This medium is intended for use as one component in a complete culture environment. EpiLife is a basal medium containing essential and non-essential amino acids, vitamins, other organic compounds, trace minerals, and inorganic salts. EpiLife does not contain antibiotics, antimycotics, hormones, growth factors, or proteins. This medium is HEPES and bicarbonate buffered and is designed for use in an incubator with an atmosphere of 5% CO₂/95% air. To support plating and long term proliferation of human keratinocytes or human corneal epithelial cells, EpiLife must be supplemented with selected hormones and growth factors. We offer a choice of growth supplements for culturing human keratinocytes and human corneal epithelial cells. Please refer to the Table on page 2 to determine the best supplement choice for your purposes. The supplements contain all of the growth factors, hormones, and extracts necessary for growth of human keratinocytes or human corneal epithelial cells in EpiLife medium.

Intended Use

EpiLife is intended for use in the routine culture of human epidermal keratinocytes or human corneal epithelial cells. Additional applications for use may include primary isolation of keratinocytes or corneal epithelial cells. Please see the table on page 2 for more information. **This product is for research use only. Not for use in animals, humans, or diagnostic procedures.**

Caution: If handled improperly, some components of this product may present a health hazard. Take appropriate precautions when handling this product, including the wearing of protective clothing and eyewear. Dispose of properly.

Storage and Stability

EpiLife is stored at 4° C in our facility and is shipped at ambient temperature. Upon receipt, the medium should be stored at 4° C and should not be frozen. **Protect from light.** Several components of this tissue culture medium are light-labile, and we recommend that the medium not be exposed to light for lengthy periods of time. If the medium is warmed prior to use, do not exceed 37° C. When stored in the dark at 4° C, the product is stable until the expiration date on the label.

Please use the supplementation instructions (page 2) to prepare the medium for use.

EpiLife is a Registered Trademark of Cascade Biologics, Inc.

EpiLife® Extended-Lifespan Culture Systems from Cascade Biologics			
	For keratinocytes:		For corneal epithelial cells:
	Undefined	Defined	Undefined
Basal Medium	EpiLife Medium (M-EPIcf-500, M-EPIcf/PRF-500, or M-EPI-500-CA)	EpiLife Medium (M-EPIcf-500, M-EPIcf/PRF-500, or M-EPI-500-CA)	EpiLife Medium (M-EPIcf-500, M-EPIcf/PRF-500, or M-EPI-500-CA)
Growth Supplement	HKGS- V2 (S-014-5)* <i>or</i> HKGS (S-001-5) <i>or</i> HKGS Kit (S-001-K)	EDGS** (S-012-5)	HCGS (S-009-5)
Coating Matrix Kit	Coating Matrix Kit (R-011-K)*	Coating Matrix Kit (R-011-K)**	N/A
Subculture Reagent	Trypsin/EDTA (R-001-100) <i>or</i> Recombinant Trypsin/EDTA (R-009-50)	Trypsin/EDTA (R-001-100) <i>or</i> Recombinant Trypsin/EDTA (R-009-50)	Trypsin/EDTA (R-001-100) <i>or</i> Recombinant Trypsin/EDTA (R-009-50)
Subculture Reagent	Trypsin Neutralizer (R-002-100)	Defined Trypsin Inhibitor (R-007-100)	Trypsin Neutralizer (R-002-100)

*For optimal cell growth, HKGS- V2 should be used in conjunction with EpiLife Medium and Coating Matrix Kit.

**For optimal cell growth, EDGS should be used in conjunction with Coating Matrix Kit and Defined Trypsin Inhibitor.

Preparation of Supplemented EpiLife Medium

Depending on your application, obtain one of the growth supplement products shown in the table above (growth supplements sold separately). To supplement one bottle of EpiLife medium with the appropriate growth supplement, we recommend the following protocol:

1. Thaw the chosen growth supplement according to the instructions provided with that product. Make sure that the caps of all of the bottles are tight.
2. Gently swirl the bottle(s) of supplement. Avoid splashing the supplement into the cap of the bottle or causing the supplement to foam.
3. Obtain one bottle of EpiLife with Calcium. Wipe the outside of the containers with a disinfecting solution such as 70% ethanol or isopropanol and move the containers to a laminar flow culture hood.
4. To add the EDGS, HKGS-V2, HKGS or HCGS, use sterile technique in a laminar flow culture hood and transfer the entire contents of the bottle of supplement to the bottle of EpiLife medium. To add the HKGS Kit, use sterile technique and transfer the desired amount of each component of the HKGS Kit to the bottle of medium.

Note: addition of less than the entire amount of any component may affect the performance of the supplemented medium. If antibiotics/antimycotics are desired, add the antibiotic/antimycotic solution included in HKGS Kit using the same technique as above.

5. Tightly cap the bottle of supplemented medium and swirl the contents to ensure a homogeneous solution. Avoid causing the medium to foam.

Storage and Stability of Supplemented EpiLife

Once EpiLife medium has been supplemented with EDGS, HKGS-V2, HKGS, HKGS Kit, or HCGS, the supplemented medium should be stored at 4° C and should not be frozen. When stored at 4° C, the supplemented medium is stable for 1 month.

Terms and Conditions of Sale

Cascade Biologics, Inc. (hereinafter, CBI) warrants that its products will perform according to the information provided in various publications that it distributes, and as described herein for the intended shelf life of the product when stored under the conditions prescribed by CBI. If any product does not perform according to the published information provided by us, CBI will replace the product free of charge to the original customer. The remedy of product replacement shall be the customers' sole and exclusive remedy for defective product, unless CBI is unable to deliver replacement product, in which case CBI shall reimburse the customer for the purchase price of the defective product. Customer understands that the foregoing limited warranty is in lieu of all other warranties and CBI hereby disclaims all other warranties including, but not limited to, implied warranty of merchantability or fitness for adequacy for any particular purpose or use. CBI shall not be liable to customer or to any party claiming through customer for any incidental or consequential damages, including, but not limited to any lost profits, lost savings, or lost business, whether arising out of contract, tort, or otherwise. Customer acknowledges that products are intended for research use only and are not to be used in animals, humans, or diagnostic procedures.