



Composition of EpiLife Defined Growth Supplement (EDGS)

Catalog # S-012-5

EDGS is a sterile 100X solution.

<u>Component</u>	<u>[EDGS]</u>	<u>[Final]</u>
Bovine Serum Albumin (>95% pure), essentially fatty acid-free (<0.005% FAs)	3.0 mg/ml	30 µg/ml
Bovine Transferrin (>95% pure)	0.5 mg/ml	5 µg/ml
rhIGF-1 (>95% pure)	1.0 µg/ml	10 ng/ml
rhEGF (>95% pure)	0.1 µg/ml	1 ng/ml
Hydrocortisone (>96% pure)	1.1 µg/ml	11 ng/ml
PGE-2 (synthetic, >98% pure)	1.8 µg/ml	18 ng/ml