



Composition of Human Corneal Growth Supplement (HCGS)

Catalog # S-009-5

HCGS is a sterile 100X solution.

<u>Component</u>	<u>[HCGS]</u>	<u>[Final]</u>
Bovine Pituitary Extract	20%	0.2%
Epidermal Growth Factor (mouse)	100 ng/ml	1 ng/ml
Hydrocortisone	18 µg/ml	0.18 µg/ml
Insulin	500 µg/ml	5 µg/ml
Transferrin	500 µg/ml	5 µg/ml