



Composition of Human Keratinocyte Growth Supplement (HKGS)
Catalog # S-001-5

HKGS is a sterile 100X solution.

<u>Component</u>	<u>[HKGS]</u>	<u>[Final]</u>
Bovine Pituitary Extract	20%	0.2%
Epidermal Growth Factor (hr*)	20 ng/ml	0.2 ng/ml
Hydrocortisone	18 µg/ml	0.18 µg/ml
Insulin	500 µg/ml	5 µg/ml
Transferrin	500 µg/ml	5 µg/ml

*hr = human recombinant