



Composition of Low Serum Growth Supplement (LSGS)

Catalog # S-003-10

LSGS is a sterile 50X solution.

<b><u>Component</u></b>	<b><u>[LSGS]</u></b>	<b><u>[Final]</u></b>
Fetal Bovine Serum	95%	1.9%
basic Fibroblast Growth Factor (hr*)	150 ng/ml	3 ng/ml
Heparin	500 µg/ml	10 µg/ml
Hydrocortisone	50 µg/ml	1 µg/ml
Epidermal Growth Factor (hr*)	500 ng/ml	10 ng/ml

\*hr = human recombinant