



SMDS

Smooth Muscle Differentiation Supplement

Catalog Number: S-008-5

5 ml

Instructions for storage and use

Product Description

Smooth Muscle Differentiation Supplement (SMDS) is a sterile, concentrated (100X) solution intended for use with basal Medium 231 for differentiation of normal human vascular smooth muscle cells. SMDS is an ionically-balanced supplement containing fetal bovine serum and heparin. Each 5 ml bottle of SMDS is the correct amount of differentiation supplement for a 500 ml bottle of basal Medium 231. When a 500 ml bottle of medium is supplemented with SMDS, the final concentrations of the components in the supplemented medium are: fetal bovine serum, 1% v/v; and heparin, 30 µg/ml.

Intended Use

SMDS is intended for use in conjunction with basal Medium 231 to induce differentiation of normal human vascular smooth muscle cells. Vascular smooth muscle cells cultured in Medium 231 supplemented with Smooth Muscle Growth Supplement (SMGS, cat. # S-007-25) will proliferate. Upon transferring these proliferating cultures into Medium 231 supplemented with SMDS, the cells will undergo a cessation of growth, a change in cellular morphology from slender stellate cells to enlarged rectangular or triangular shaped cells, and will express increased amounts of smooth muscle α -actin over the course of 2 to 8 days.

Reducing serum concentration and adding heparin to smooth muscle cell cultures are both commonly used procedures in differentiation studies of smooth muscle cells *in vitro*. For example see:

Newcomb, P.M. and Herman, I.M.; Pericyte growth and contractile phenotype: Modulation by endothelial-synthesized matrix and comparison with aortic smooth muscle. *J. Cell Physiol* 155: 385-393; 1993

Reilly, C.F., Fritze, L.M.S., Rosenberg, R.D.; Heparin inhibition of smooth muscle proliferation: a cellular site of action. *J. Cell Physiol* 129: 11-19; 1986

Castellot, J.J., Favreau, L.V., Karnovsky, M.J., Rosenberg, R.D.; Inhibition of vascular smooth muscle cell growth by endothelial cell-derived heparin. *J. Biol Chem* 257(19): 11256-11260; 1982

This product is for research use only. Not for use in animals, humans, or diagnostic procedures.

Caution: If handled improperly, some components of this product may present a health hazard. Take appropriate precautions when handling this product, including the wearing of protective clothing and eyewear. Dispose of properly.

Storage and Stability

SMDS is stored at -20° C at our facility and is shipped on dry ice. Upon receipt, the product should be stored at -20° C in a freezer that is not self-defrosting. When stored at -20° C the product is stable until the expiration date shown on the label.

Thawing

To thaw, place the product in a 37° C water bath or overnight at 4° C. If thawed in a water bath, do not leave the product at 37° C after the product has thawed. For instructions on adding SMDS to basal Medium 231, please refer to the instructions that accompany the basal medium.

Terms and Conditions of Sale

Cascade Biologics, Inc. (hereinafter, CBI) warrants that its products will perform according to the information provided in various publications that it distributes, and as described herein for the intended shelf life of the product when stored under the conditions prescribed by CBI. If any product does not perform according to the published information provided by us, CBI will replace the product free of charge to the original customer. The remedy of product replacement shall be the customers' sole and exclusive remedy for defective product, unless CBI is unable to deliver replacement product, in which case CBI shall reimburse the customer for the purchase price of the defective product. Customer understands that the foregoing limited warranty is in lieu of all other warranties and CBI hereby disclaims all other warranties including, but not limited to, implied warranty of merchantability or fitness for adequacy for any particular purpose or use. CBI shall not be liable to customer or to any party claiming through customer for any incidental or consequential damages, including, but not limited to any lost profits, lost savings, or lost business, whether arising out of contract, tort, or otherwise. Customer acknowledges that products are intended for research use only and are not to be used in animals, humans, or diagnostic procedures.